



HAMILTON-MR1

Intelligent Ventilation from the ICU to MRI





Intelligent Ventilation from the ICU to MRI

HAMILTON-MR1 – The MR Conditional ventilator

The fully featured ICU ventilator, HAMILTON-MR1, guarantees uncompromised continuous ventilation care from the ICU to the MRI scanner and back. Its reliability and high performance, with advanced lung-protective strategies and patient-adaptive modes, make the HAMILTON-MR1 the ideal choice for any critical care department that needs to transport ventilated patients to the MRI department.

- ✓ MR Conditional (up to 50 mT)
- ✓ Integrated TeslaSpy gaussmeter
- ✓ Adult, pediatric, and neonatal ventilation
- ✓ More than 9 hours of battery operating time
- ✓ Independence from gas cylinders or compressors
- ✓ Noninvasive ventilation and integrated high flow oxygen therapy*
- ✓ Advanced ventilation modes including ASV® - Adaptive Support Ventilation



You can take your patients from the ICU down to the MRI suite and not have to change a thing about the ventilation, even when they are on an advanced mode. That is a true advantage because you are not risking lung derecruitment and a patient setback, which would keep him in the hospital longer and make it more uncomfortable for him.

Dr. Thomas Berlin, Director of Respiratory Care
Florida Hospital Orlando, Orlando (FL), USA

Adapted to the MR environment

Integrated gaussmeter for more safety


Positioning a medical device too close to the MRI scanner can have fatal consequences. The integrated gaussmeter, TeslaSpy, continuously monitors the magnetic field and gives you an audible and visual signal if you are getting too close. For maximum safety, TeslaSpy continues monitoring even when the ventilator is switched off.

Close to the patient

The HAMILTON-MR1 is the first ventilator able to be used at a magnetic field strength of 50 mT (equivalent to 1 m distance for a 3T static magnetic field scanner), without creating any MR image artifacts.

Ideal for clinical transport

The small and rugged housing of the HAMILTON-MR1 makes it easy to handle and optimal for clinical transport. With its integrated turbine, full range of modes, and powerful internal batteries, the HAMILTON-MR1 accompanies your patient from the ICU to MRI and back, providing uncompromised ventilation in a compact design.



TeslaSpy continuously monitors the magnetic field and gives you an audible and visual signal if you are getting too close



Intrahospital transport to the MRI

Optimal performance

The right ventilation mode for your patients

In addition to conventional and modern modes of invasive and noninvasive ventilation, the HAMILTON-MR1 also offers the option of an integrated high flow oxygen therapy mode. Using the same device and breathing circuit, you can change the interface in just a few quick steps and adjust the therapy to best meet your patients' needs. This guarantees that before, during, and after the MRI procedure, your patients receive the same high level of ventilation care as at the bedside.

Adaptive synchronization

The IntelliTrig function automatically adjusts the inspiratory and expiratory trigger sensitivity to potential leaks and ensures adaptive synchronization with the patient's breathing pattern. This is achieved both for invasively and noninvasively ventilated patients.

Give your patients a voice

In pressure-controlled modes (PCV+, SPONT, PSIMV+), an optional feature enables use of conventional speaking valves with the HAMILTON-C1. Monitoring, triggering, and alarm management have been adjusted to allow the use of speaking valves.



Invasive ventilation





High flow oxygen therapy

Product overview

- 1 Integrated TeslaSpy
- 2 Patient interfaces and ports
- 3 Press-and-turn knob
- 4 Ventilation Cockpit
- 5 360° visible alarm lamp
- 6 Power supply and connectors
- 7 Optional handle







Patients in the medical intensive care unit could be extubated earlier following the introduction of ASV.⁵

Time to Extubation Readiness (d)

More safety and comfort for your patients

Enhanced patient comfort

Every Hamilton Medical ventilator features the intelligent ventilation mode ASV (Adaptive Support Ventilation). ASV measures the patient's lung mechanics and activity on a breath-by-breath basis and automatically adjusts ventilation, from intubation to extubation. Since its introduction in 1997, ASV has become well established in intensive care units and has been shown to improve patient/ventilator interaction.^{1), 2)}

Lung-protective ventilation

ASV ensures via an optimal breathing pattern that the patient receives the set minute volume, irrespective of the patient's activity. As part of this process, ASV employs lung-protective strategies to minimize complications from AutoPEEP and volutrauma/barotrauma. ASV also prevents apnea, tachypnea, excessive dead-space ventilation, and excessively large breaths.³⁾

Decreased ventilation time

Clinical studies show that:

- ASV supports the earliest possible spontaneous breathing by the patient^{4), 5)}
- ASV shortens the ventilation time in various patient groups^{4), 5)}

1 Iotti GA. Intensive Care Med. 2010 Aug;36(8):1371-9. | 2 Sulzer CF. Anesthesiology. 2001 Dec;95(6):1339-45.
3 Sulemanji D. Anesthesiology. 2009 Oct;111(4):863-70. | 4 Kirakli C. Eur Respir J. 2011 Oct;38(4):774-80 | 5 Chen CW. Respir Care. 2011 Jul;56(7):976-83.

Ease of use

Intuitive operation


In close cooperation with users and ventilation experts, our engineers have designed the HAMILTON-MR1 user interface to allow intuitive operation and direct access to important settings. All Hamilton Medical ventilators are operated according to the same principles, which makes switching between different devices very easy.

Easy-to-understand monitoring

Ventilators display large amounts of data that is often difficult to interpret. The configurable touch screen display, referred to as the Ventilation Cockpit, consolidates the diverse monitoring data, and presents it numerically and in various graphics panels. These easy-to-understand views provide an at-a-glance overview of the patient's current ventilation status, and offer a reliable basis for therapy decisions.

More time for your patients


In ASV mode, the ventilator continuously adjusts to the patient's breathing activity and lung conditions. This means fewer user interactions are required and fewer alarms are generated¹⁾, giving you more time for your patients.



Monitoring window






More time for your patient



The Ventilation Cockpit

- ① Dynamic Lung - Provides a real-time display of lung compliance, resistance, breathing activity, and pulse rate
- ② Direct access to the most important settings
- ③ The four most important monitoring parameters
- ④ Configurable waveforms for flow and pressure
- ⑤ Display options of the Ventilation Cockpit:
 - a) ASV Graph
 - b) Vent Status
 - c) Trends (not shown)
 - d) Loops (not shown)





“

Until now, we ventilated our ICU patients during MRI scanning with an anesthesia ventilator. Therefore, we had to consult an anesthetist every time to handle the equipment. With the HAMILTON-MR1, we are now completely independent.

Dr. Adrian Wäckerlin, Head of Intensive Care
Cantonal Hospital Grisons, Chur, Switzerland

Increased efficiency

Integrated commercial considerations

Ventilators are capital goods that need to be evaluated for cost efficiency. Factors including treatment costs and the use of human resources play an important role in this process. Assembled with an extensive standard equipment package that is easy to maintain, Hamilton Medical ventilators are an attractive investment with respect to purchase price and operating costs.

Reduction of treatment costs

For each day where ventilation is no longer required, treatment costs are reduced by 1,500 USD on average.¹⁾ It has been shown that the use of Hamilton Medical ventilators and ASV can reduce ventilation time. In addition, the ventilator is then available for the next patient much earlier. A shorter ventilation time also reduces the risk of ventilator-associated pneumonia (VAP), which can result in costs of up to 57,000 USD per case.²⁾

Better use of human resources

Hamilton Medical ventilators, along with ASV, can reduce the time needed for standard settings and alarm management while maintaining ventilation quality.^{3), 4)} This frees up time for other aspects of patient care. Thanks to the ease of operation, consistent operating concepts across devices, and the free e-learning offerings from Hamilton Medical, the effort for education and training is also reduced.

1 Dasta JF et al. Critical Care Med. 2005 Jun;33:1266-71 | 2 Cocanour CS et al. Surg Infect. 2005 Spring;6:65-72
3 Iotti GA. Intensive Care Med. 2010 Aug;36(8):1371-9 | 4 Petter AH. Anesth Analg. 2003 Dec;97(6):1743-50

Attention to detail

Useful features to make your life easier

The trolley's auto-lock brake locks the wheels as soon as you release the handle to prevent it rolling accidentally towards the MR scanner. Hooks on each side of the trolley enable convenient storage of the breathing circuit and oxygen hose. The optional kit with quick-lock function and specially-designed handle allows you to remove the device at the press of a button and attach it directly to the bed.

Customizable user interface

You can configure the display layout with different waveforms, loops, trends, or intelligent panel graphics to suit your institution's needs and protocols. Nurses and clinicians can have their own preferred layout. Access the Monitoring window with the touch of a button at any time during active ventilation.

Free and open e-learning on mechanical ventilation

Join over 22,000 users on the Hamilton Medical College e-learning platform. It provides free and open e-learning modules on the basics of mechanical ventilation, as well as on Hamilton Medical products and features. Register now at college.hamilton-medical.com.

For some modules, a certificate is issued upon successful completion. You can even receive Continuing Respiratory Care Education (CRCE) credits from the American Association of Respiratory Care (AARC) for some modules.



Auto-lock brake



Clinicians using the e-Learning platform

Neonatal ventilation

Tidal volumes as low as 2 ml

With the neonatal option, the HAMILTON-MR1 provides tidal volumes as low as 2 ml for effective, safe, and lung-protective ventilation even for the smallest patients.¹⁾ The proximal flow sensor specifically developed for neonates precisely measures the pressure, volume, and flow directly at the infant's airway opening, ensuring the required trigger sensitivity. This provides improved synchronization and less work of breathing.

Adaptive synchronization, even with uncuffed tubes

Leaks are one of the issues encountered in the ventilation of neonates, as a result of using uncuffed tubes. The IntelliTrig leak compensation function automatically adjusts the inspiratory and expiratory trigger sensitivity to potential leaks. This enables adaptive synchronization with the neonate's breathing pattern.

nCPAP - Automatic adaptation, fewer interventions

The HAMILTON-MR1's nCPAP mode is designed in such a way that you only need to set the desired CPAP pressure. The flow is subsequently adjusted automatically based on the patient condition and potential leaks. This prevents unintended peak pressures and guarantees highly efficient leak compensation. Flow adjustment occurs very rapidly due to near-patient pressure measurement and the high sensitivity of the measurement.



Effective, safe, and lung-protective ventilation for the most fragile patients



Neonate with nCPAP mask

Hamilton Medical

Intelligent Ventilation since 1983

In 1983 Hamilton Medical was founded with a vision: To develop intelligent ventilation solutions that make life easier for patients in critical care and for the people who care for them. Today, Hamilton Medical is a leading manufacturer of critical care ventilation solutions for a wide variety of patient populations, applications, and environments.

The right ventilation solution for every situation

The ventilators from Hamilton Medical ventilate all of your patients; in the intensive care unit, during an MRI procedure and in all transport situations, from the neonate to the adult. Each of these ventilators is equipped with the same standardized user interface and uses the same Intelligent Ventilation technologies. This enables Hamilton Medical ventilators to help you to:

- ✓ Increase the comfort and safety of your patients
- ✓ Make life easier for the caregivers
- ✓ Increase efficiency and return on investment

