

MEDIVATORS™  
**ENDO GATOR™**  
Irrigation Pump



# Automate Irrigation and Reduce Procedure Time

MEDIVATORS™ ENDOGATOR™ Irrigation System eliminates manual syringing. Irrigation is provided via footswitch, reducing procedural time and increasing visibility.

## Pump Specifications

### Electrical

Voltage/Current Requirement:  
120V/240V~50-60Hz

### Pump

Flow Rate: 0-550 ml/min  
Pressure: 0-48 psi

### Ideal Flow Settings

Auxiliary Water Channel: 40-50  
Biopsy Channel: 70 to Max

### Physical Dimensions

Height: 152 mm    Depth: 203 mm  
Width: 178 mm    Weight: 2.3 kg

## MEDIVATORS™ ENDOGATOR™ Accessories

### Safety












- Reduced risk of cross contamination
- Endoscope connectors also available in single patient version

### Compatibility

- Works with all sterile water bottles
- Connects to all major brands of endoscopes

### Efficiency

- Minimizes reprocessing requirements
- Creates consistent practices
- Cost effective and completely disposable
- Endoscope connectors for all methods of irrigation (Biopsy, Jet)

	Tubing	Biopsy Channel Irrigation	Jet Channel Irrigation	
<b>OLYMPUS®</b> GI Endoscopes	200230U 	100303 	100115 	100241 
<b>PENTAX®</b> GI Endoscopes		100136 	100116 	100242 
<b>STORZ FUJINON®</b> GI Endoscopes	200230FU 	100303 	Backflow Valve 100126 	Direct Luer Connection on Tubing 

ENDO GATOR™ and MEDIVATORS™ are trademarks of Medivators Inc.  
PENTAX® is a registered trademark of Hoya Corporation.  
OLYMPUS® is a registered trademark of OLYMPUS Corporation.  
FUJINON® is a registered trademark of FUJIFILM Corporation.



MEDIVATORS Inc. Headquarters  
14605 28th Avenue North  
Minneapolis, MN  
55447-4822 USA  
Tel: +1.800.444.4729  
Fax: +1.763.553.3387

MEDIVATORS BV  
Sourethweg 11  
6422 PC Heerlen  
The Netherlands  
Tel: +31.45.5.471.471  
Fax: +31.45.5.429.695

MEDIVATORS Asia/Pacific Pte. Ltd.  
1 A International Business Park  
#05-01 Singapore 609933  
Tel: +65.6227.9698  
Fax: +65.6225.6848

MEDIVATORS Beijing Representative Office  
Room 708, 7th Floor Kaiheng Center, Block B  
No. 2 Chaoyangmennei Street  
Dongcheng District, Beijing China 100010  
Tel: +8610.6567.8446  
Fax: +8610.6567.8445

[www.medivators.com](http://www.medivators.com)

© MEDIVATORS Inc. P/N 50098-506