


BRAVOS inspires new victories in the fight against cancer



It's a brave new world
in brachytherapy





BRAVOS lets you focus on your patients

When you have confidence in your tools, it instills confidence in your patients. That's why Varian has combined more than 50 years of our experience in brachytherapy with your clinical experience to design an intuitive afterloader that enables you to focus on your patients.

The BRAVOS advantage

Pre-treatment checklist

Pre-treatment checklist designed for better patient safety and treatment efficiency, including e-Signature prior to therapy.

Patient safety

Time-out for patient safety requires a comprehensive review before a procedure to ensure all protocols are accounted for.


Simplified use

LED display confirms applicator connections with easy-to-read colored indicator lights. Transfer guide tubes can be the same length and provide simple connection to be used in any channel.

Integration

To track a patient's progress from diagnosis to post-treatment review, BRAVOS integrates with the BrachyVision™ 3D treatment planning system, the ARIA® oncology information system (OIS), and the Eclipse™ treatment planning system.

The BRAVOS™ afterloader system is the result of clinically driven design, developed with the help of professionals like you—radiation oncologists, medical physicists, therapists, dosimetrists, service engineers, and IT specialists. With your input, we designed and engineered BRAVOS to be efficient, integrated, and easy to use. Welcome to the brave new world of BRAVOS.




Efficient

- User interface streamlines the overall treatment from planning to delivery.
- Customizable pre-treatment checklist, similar to surgical procedures, can be set to require checklist review and an e-Signature before proceeding, instilling greater confidence and accountability.
- Time-out for patient safety requires you to stop and review all brachytherapy treatment protocols before proceeding.
- Recalibration of source position by the clinic can optimize uptime. The BRAVOS system with the CamScale™ automated position verification system provides tools to visualize and recalibrate both source and dummy positions between source exchanges, minimizing the need for a service call.
- Automatic length measurement allows accuracy between treatment planning and treatment delivery, regardless of the treatment length.

Integrated

- Integrates with the BrachyVision 3D treatment planning system so you can create comprehensive treatment plans for use with BRAVOS.
- Varian OIS interoperability means BRAVOS works seamlessly with ARIA and Eclipse to manage the patient's entire journey—from initial diagnosis through post-treatment follow-up.
- Connects with the ARIA oncology information system to effectively manage the administrative aspects of your department.
- Compatible with the Varian family of applicators.



Easy to use

- Touchscreen controls provide fast, easy access to patient treatment information in the vault.
 - The SmartConnect® remote monitoring service allows BRAVOS systems to be monitored by Varian's service engineers who can troubleshoot and easily resolve issues in real time.
 - Predictive maintenance enables you to proactively maximize your uptime.
- * Fighting cancer takes more than a machine. It's us working together. With your help, we've created a connected ecosystem designed around you, your clinic, and your patients. Working together, we can find new and unexpected ways of elevating treatment and delivery to reach your next victories. Let's inspire a brave new world of brachytherapy, together. Come join us.

Brachytherapy, an integral part of the evolving Varian ecosystem.

Learn more at www.varian.com/brachytherapy



Not all products or features are available in all markets.

Important Safety Information

Radiation therapy is not appropriate for all cancers. For more information, please visit www.varian.com/safety.

Medical Advice Disclaimer

Varian as a medical device manufacturer cannot and does not recommend specific treatment approaches. Individual treatment results may vary.

varian

varian.com

**USA, Corporate
Headquarters and
Manufacturer**

Varian Medical Systems, Inc.
3100 Hansen Way
Palo Alto, CA 94304
Tel: 650.424.5700
800.544.4636

**Headquarters Europe,
Eastern Europe, Middle &
Near East, India, Africa**

Varian Medical Systems
International AG
Steinhausen, Switzerland
Tel: 41.41.749.8844

Asia Pacific Headquarters

Varian Medical Systems
Pacific, Inc.
Kowloon, Hong Kong
Tel: 852.2724.2836

Australasian Headquarters

Varian Medical Systems
Australasia Pty Ltd.
Sydney, Australia
Tel: 61.2.9485.0111

Latin American Headquarters

Varian Medical Systems
Brasil Ltda.
São Paulo, Brazil
Tel: 55.11.3457.2655

© 2019, 2020 Varian Medical Systems, Inc. Varian, Varian Medical Systems, ARIA and SmartConnect are registered trademarks, and BrachyVision, BRAVOS, CamScale and Eclipse are trademarks of Varian Medical Systems, Inc.